

Mycologic

User Guide

<https://web-production-77f82.up.railway.app/>

Quick Check

A fast, one-off health reading from a single photo. No account needed.

1. Open the Quick Check tab
2. Upload a photo of your specimen
3. Click "Run analysis"
4. View the annotated result and health score (0–100)

The score is based on the image only — green = healthy, yellow = dry/aging, red = contamination, grey = exposed substrate.

Sample Archive

For long-term tracking. Create a profile for each specimen and log observations over time. Requires a free account.

1. Go to "My Samples" → click "+ New sample"
2. Give it a name and optional cover photo
3. Inside the sample, upload detection photos — each one is scored and added to the timeline
4. Upload sensing data (CSV from ESP32/BME280) to add humidity into the score
5. Upload CloudCompare scans and humidity maps to build a complete record

The score here combines three dimensions: Image (60 pts) + Humidity (30 pts) + Trend (10 pts). Below 60 = the system flags that intervention may be needed.

Contribute

Help improve the detection model by contributing annotated images. No account needed.

1. Go to the "Contribute" tab
2. Upload a photo of any mycelium surface
3. Draw polygons around features and label them (healthy, dry aging, contamination, or exposed substrate)
4. Click "Submit" — your annotation is added to the shared training dataset

Every contribution helps the model learn. The more diverse the data, the better the detection.